



dietary
supplements



biocidal
products



medicinal
products



semen
extenders



veterinary
products

BIOSOLWENS LUX[®] - boar semen extender

It enables to store diluted boar semen **up to 5-8 days** at room temperature (15-20°C).

CHARACTERISTICS:

- acrosomal integrity
- plasma membrane functionality
- higher motility
- higher conception rate

COMPONENTS:

The extender contains carefully chosen components, which preserve and protect the acrosome and cell membrane of spermatozoa against damage. It contains antibiotic – lincomycin.

Its properties and efficiency have been tested under laboratory and farm conditions by cooperation with The Polish National Research Institute of Animal Production in Balice (Kraków).

PREPARATION:

Biosolwens extender is simple to use. You have to dissolve the content of one bag in an adequate amount of distilled water at a temperature of +30°C to +40°C (stir gently). Diluted semen should be store and transport at temperature from +15°C to +20°C.

STORAGE:

Unopened bags should be stored in a dry and dark place.

FORMATS

Sachets with powder to dissolve in: 1 liter, 5 l, 10 l, 25 l, 50 l, or 100 l

EXPIRATION:

1 year from packaging date stated on the sachet.

MANUFACTURER:

Pharmaceutical Research and Production Plant “BIOCHEFFA” Sp. Z o.o.

ul. Kasztanowa 3, 41-205 Sosnowiec, POLAND

tel/fax: +48 32 2916968

person to contact: Artur Biliński

e-mail: marketing@biocheffa.pl

www.biosolwens.pl

www.biocheffa.pl

