



dietary  
supplements



biocidal  
products



medicinal  
products



semen  
extenders



veterinary  
products

## BIOSOLWENS PLUS<sup>®</sup> - boar semen extender

It enables to store diluted boar semen **up to 7-10 days** at room temperature (15-20°C).



### CHARACTERISTICS:

- acrosomal integrity
- plasma membrane functionality
- higher motility
- higher conception rate

### COMPONENTS:

The extender contains carefully chosen components, which preserve and protect the acrosome and cell membrane of spermatozoa against damage.

It contains antibiotic – **gentamicin**.

Its properties and efficiency have been tested under laboratory and farm conditions by cooperation with The Polish National Research Institute of Animal Production in Balice (Kraków).

### PREPARATION:

Biosolwens extender is simple to use. You have to dissolve the content of one bag in an adequate amount of distilled water at a temperature of +30°C to +40°C (stir gently). Diluted semen should be store and transport at temperature from +15°C to +20°C.

### STORAGE:

Unopened bags should be stored in a dry and dark place.

### FORMATS:

Sachets with powder to dissolve in: 1 liter, 5 l, 10 l, 25 l, 50 l, or 100 l

### EXPIRATION:

1 year from packaging date stated on the sachet.

### MANUFACTURER:



Pharmaceutical Research and Production Plant “BIOCHEFFA” Sp. Z o.o.

ul. Kasztanowa 3, 41-205 Sosnowiec, POLAND

tel/fax: +48 32 2916968

person to contact: Artur Biliński

e-mail: [marketing@biocheffa.pl](mailto:marketing@biocheffa.pl)

[www.biosolwens.pl](http://www.biosolwens.pl)

[www.biocheffa.pl](http://www.biocheffa.pl)